


Prodigy ITFAP

Quick Guide to Identification of Thin Films From Angular Profiles



SPECS Surface Nano Analysis GmbH

Voltastrasse 5

13355 Berlin

Germany

Tel: +49 30 46 78 24-0

Fax: +49 30 46 42 08-3

Email: support@specs.com

Web: www.specs-group.com



SPECS User Manual

Prodigy ITFAP—Quick Guide to Identification of Thin Films From Angular Profiles

SpecsLab Prodigy Version 4.103.1

June 6, 2023

SPECS order number reference: N/A

©2023. All rights reserved. No part of this manual may be reproduced without the written permission of SPECS Surface Nano Analysis GmbH. Disclosing this instruction to third parties is only permitted with the prior consent of SPECS Surface Nano Analysis GmbH. Other product and company names mentioned in this document may be the trademarks or registered trademarks of their respective owners and are used for identification purposes only.

Table of Contents

Chapter 1 – Overview

1.1	User Interface	1
1.2	Key Features	2

Chapter 2 – Prerequisites

2.1	General Information	3
2.2	Configuration of System Geometry	3
2.2.1	Orientation of the Analyzer	5
2.2.2	Orientation of the X-ray Source	6
2.3	Angle-resolved Transmission Function	7

Chapter 3 – ARXPS Depth Profile Analysis

3.1	Measurement of Detailed Angle-resolved XPS Spectra	9
3.2	Peak Fitting	9
3.3	Depth Profile Analysis	10
3.4	Sample Model Definition	12
3.5	Sample Model Presets	18


This page intentionally left blank.

Chapter 1 – Overview

ITFAP is a software module integrated in SpecsLab Prodigy for non-destructive depth profile analysis of layered thin film samples by ARXPS.

1.1 User Interface

In SpecsLab Prodigy, the Depth Profile Analysis module appears with the interface shown below:



1.2 Key Features

The ITFAP module is a separately licensed option integrated with SpecsLab Prodigy. It compares angle profiles—peak areas as a function of emission angle—measured by parallel ARXPS and calculated for a layer model of the sample.

The following features are provided with ITFAP:

- Seamless integration with data acquisition
- Immediate layer structure information from relative depths of atomic species
- Easy input of proposed sample layer model
- Instantaneous optimization of sample model layer thicknesses
- Evaluation of consistency between sample model and ARXPS measurement via goodness of fit measure
- Quick comparison of sample models via sample model preset system
- Extensive parameter sourcing and control via SpecsLab Prodigy Chemical Databases

Chapter 2 – Prerequisites

2.1 General Information

The ITFAP module is recommended to users who have suitable experience in the field of ARXPS.

For advice and assistance, please contact SPECS support:

Tel. +49 30 46 78 24-0

email: support@specs.com

2.2 Configuration of System Geometry

Non-destructive depth profile analysis by ARXPS requires correct configuration of the system geometry. The SpecsLab Prodigy Configuration Tool is used to specify the orientation of

- analyzer, and
- X-ray source.

Both orientations are specified within a coordinate system attached to the analysis chamber (system coordinate system). For depth profile analysis, the orientation of the sample surface normal within the system coordinate system must also be specified. This is done during data analysis in ITFAP itself, not within the configuration tool.


Note

It is possible to freely choose the orientation of the system coordinate system relative to the analysis chamber. Once this choice is made, it needs to be maintained when specifying the analyzer, source, and sample orientation.


Often, the system coordinate system is chosen such that the sample surface plane lies within the x-y plane for a typical sample orientation.

To specify analyzer and source orientation, open the respective device configuration in the configuration tool and complete the steps described in the following paragraphs.

The figure below shows the analyzer device configuration:




To access the source device configuration, switch from the analyzer to the X-ray source:



2.2.1 Orientation of the Analyzer

The below screen of the configuration tool shows the angular parameters (1) specifying the orientation of the analyzer with respect to the system coordinate system:

- *Polar Angle*
- *Azimuth Angle*
- *Rotation Angle*



The geometry sketch (2) in the image describes these parameters.

Coordinate Systems

In the above figure, the system coordinate system is colored black. It consists of the axes x, y, and z. The analyzer coordinate system consists of the axes

- X_{Analyzer}
- Y_{Analyzer}
- Z_{Analyzer}

The analyzer coordinate system is a right-handed coordinate system, where Z_{Analyzer} is the optical axis of the analyzer lens and X_{Analyzer} points in the direction of increasing electron kinetic energy on the detector.

How to Specify the Analyzer Orientation

Starting from a situation where the analyzer and system coordinate system are congruent, the analyzer orientation results from the following successive right-handed rotations:

1. Rotation of the analyzer by the *Polar angle* around the y-axis.
2. Rotation of the analyzer by the *Azimuth angle* around the z-axis.
3. Rotation of the analyzer by the *Rotation angle* around the z_{Analyzer} -axis.

2.2.2 Orientation of the X-ray Source

In the Configuration Tool, the source orientation is specified with respect to the system coordinate system using the angular parameters:

- *Polar Angle*
- *Azimuth Angle*

☒

Unique Device Name: X-Ray Dummy

Device Type: XRayDummy

Visible Name X-Ray Dummy

Device Settings

Anode 1 Al

Anode 2 Mg

Geometry (ISO)

Polar Angle 54.50 °

Azimuth Angle 90.00 °

Communication

Type Dummy

Revert

Delete

Further information:


Geometry:

The azimuth angle is measured counter-clockwise from the x-axis.

These parameters define the direction of the beam axis within the system coordinate system.

2.3 Angle-resolved Transmission Function

It is recommended to measure an angle-resolved transmission function specific to the geometry used for ARXPS data acquisition:



In SpecsLab Prodigy's *Transmission Correction* view, the orientation of the sample surface normal is specified via its polar and azimuth angle with respect to the system coordinate system. The analyzer and source orientation are used as specified in the configuration tool.

This page intentionally left blank.

Non-destructive depth profile analysis is performed in the following steps:


- These steps are described in more detail below.

In the first step, acquire detailed ARXPS spectra with the peaks of the atomic species to be considered for depth profile analysis. Use an angle-resolved transmission function for the given sample position and orientation.

Fit the peaks, separating the angular channels of the spectrometer.


Prodigy ITFAP 4.103.1 | June 6, 2023

SpecsLab Prodigy's *Peak Fit* operation automatically propagates peak function parameters starting from a user-selectable start channel.




3.3 Depth Profile Analysis

Apply a *Depth Profile Analysis* operation on the spectrum group.



The angular profiles and relative depths are available immediately in *Depth Profile Analysis*:


11



3.4 Sample Model Definition


The next step is to define the sample model based on the relative depths and known sample properties. For this purpose, the sample model editor is used.

To specify the substrate composition, open [Component Settings](#) by clicking here:




Enter a molecular formula (1) and assign peaks to the component (2):


13



Add layers above the substrate ...




... and specify their components:




For the sample model given, angle profiles are calculated immediately. The layer thickness is optimized for best fit of the calculated angle profiles to the measured angle profiles.


Angle profiles are only calculated for peaks assigned to components of the sample model. A comparison of *relative* calculated and measured angle profiles requires only those measured profiles to be considered, for which also a calculated profile exists. Check [Hide Unassigned Profiles](#) to exclude measured profiles without a calculated counterpart.



In the example used, relatively high uncertainty of the peak areas is apparent at the largest emission angles. These peak areas can be excluded from the fit by using the upper *Channel Selection* limit:




In the way described, more layers can be added in order to model the expected sample structure. Upon definition of a new layer, layer thicknesses are instantly optimized for the layer sequence defined. A reduced chi-squared (Red. χ^2) is calculated as a measure of goodness of fit:



3.5 Sample Model Presets

Sample models can be stored as presets:



The screenshot displays the 'Depth Profile Analysis' window in the SPECS software. The window is divided into several sections:

- System Geometry:** Contains settings for Sample, Analyzer, and Source, each with Polar and Azimuth angles. The Sample Polar angle is 0.00°, Azimuth is 0.00°. The Analyzer Polar angle is 50.00°, Azimuth is 90.00°. The Source Polar angle is 35.00°, Azimuth is 180.00°. The Rotation is -90.00°.
- Channel Selection:** Shows Lower and Upper channel ranges. Lower is 0, Upper is 16.
- Buttons:** Includes 'Hide Unassigned Profiles' and a dropdown menu set to '<none>'. A red box highlights a button in the bottom right corner of the window.

Below the 'Depth Profile Analysis' window, the 'Save As' dialog box is open. It shows the file path 'This PC > Documents > SpecsLab Prodigy > ITFAP'. The file name is 'SiO2_HfO2_SiO2_on_Si.sls' and the save type is 'SpecsLab Sample Model (*.slsm)'. The 'Save' button is highlighted.


Below the 'Save As' dialog box, the 'Depth Profile Analysis' window is shown again, but now displaying a plot and a table of results.

Plot: The plot shows 'Relative Normalized Peak Area [%]' on the y-axis (0 to 60) versus 'Emission Angle [deg]' on the x-axis (20.00 to 70.00). The plot includes data points for SiO2 (red circles), HfO2 (green circles), and Si (blue circles). A red line represents the fit for SiO2, and a green line represents the fit for HfO2. The plot also shows a legend for the data series.

Table: The table lists the layers and their thicknesses in nanometers (nm). The layers are SiO2 (0.24 nm), HfO2 (0.63 nm), SiO2 (2.71 nm), and Si (1.61 nm). The bottom layer is labeled 'Substrate'.

Legend: The legend at the bottom of the window identifies the data series: C1s (red), Hf4d3/2 (oxide) (green), Hf4d5/2 (oxide) (blue), Hf4f5/2 (oxide) (cyan), Hf4f7/2 (oxide) (magenta), Si2p (bulk) (pink), and Si2p (oxide) (purple).

Presets make comparison of different sample models particularly easy. In the screenshot below, a preset was loaded, where the sample model does not comprise an oxide layer on the silicon substrate. The larger resulting chi-squared shows immediately, that omitting the oxide layer reduces fit quality:



This page intentionally left blank.